

Diagnostic Ultrasound System

DC-6Vet



DC-6Vet is dedicated to provide veterinary ultrasound solutions meeting your daily clinical needs. DC-6Vet produces high-sensitivity color Doppler performance and high-resolution 2D images. Various functions such as THI, FTO, iTouch and iStation make your ultrasound diagnosis much more reliable. The system's specific veterinary application modules include reproductive software packages, animal-use transducer series, user-friendly operation interface and so on. All DC-6Vet outstanding solutions make it an ideal choice for any veterinary sonographer.

DC-6Vet

Diagnostic Ultrasound System



- 2D, M, PW, HPRF, Color, Power and DirPower
- Abundant software packages, including specially designed reproductive software for Dog, Cat, Equine, Bovine and Ovine
- Wide band multi-frequency transducers ranging from 2.0MHz to 12.0MHz
- Specific high frequency micro-convex veterinary transducer 6C2 (5.0/6.5/8.0 MHz) for small animal use
- THI (Tissue Harmonic Imaging) for technically difficult scanning
- FTO (Fine Tissue Optimization)
- iTouch™: Intelligent image optimization with only one button
- iStation™: Intelligent patient management platform greatly increasing your work efficiency
- Outstanding solutions for information management: CD-RW, four USB ports and optional DICOM3.0



Micro-convex array transducer
6C2 (5.0/6.5/8.0MHz)



Convex array transducer
3C5A (2.5/3.5/5.0/H4.6/H6.0MHz)



Linear array transducer
7L4A (5.0/7.5/10.0MHz)



Micro-convex array transducer
3C1 (2.5/3.5/5.0/H4.6/H6.0MHz)



Linear array transducer
10L4 (8.0/10.0/12.0MHz)



Linear array transducer
7L6 (5.0/7.5/10.0MHz)



Micro-convex endocavity transducer
6CV1 (5.0/6.5/8.0MHz)



Phased array transducer
2P2 (2.0/2.5/3.0/H3.5/H4.0MHz)



Distributed by:
Smart Works, Inc
4164 Austin Bluffs Pkwy. #361
Colorado Springs, CO 80918
1-866-999-2656

<http://www.smartworksimaging.com>
sales@smartworksimaging.com

